


# STREAMWISE




## ProCap newsletter

Version 3.2 release

### ProCap 3.2 release


New visualization: Project measurement data on model surfaces  
Introducing the new, interactive user interface



### Surface flow maps

The new projection feature maps measured data (e.g. velocity or static pressure) onto the model surface in real time and greatly enhances analysis capabilities.


[Professional only]



### Fast, intuitive interactive setup

All new interactive user interface with immediate preview and reversible changes for fast and intuitive setups. Developed based on customer feedback!

[Compact and Professional]



### Improved visualization

Improved data analysis and visuals with adjustable transparency of isosurfaces (Professional only) and model geometries.

[Compact and Professional]

# New in ProCap 3.2



## Professional

- ✓ **New feature: Surface projections**  
Real-time mapping of measured data on the model surface
- ✓ **Interactive user interface**  
Intuitive & interactive setup, direct preview and reversible changes
- ✓ **Improved visualization**  
Transparent model geometries, isosurfaces and line feature for streamlines
- ✓ **Simplified model positioning**  
Relative positioning in model coordinates
- ✓ **Efficient workflow**  
Clone the last setup as starting point for a new measurement
- ✓ **Probe & data handling**  
Listing of active COM-ports, user-defined threshold on incoming data and addition of data to stored files



## Compact

- ✓ **Interactive user interface**  
Intuitive interactive setup, direct preview and reversible changes
- ✓ **Improved visualization**  
Transparent model geometries and scalar planes
- ✓ **Simplified model positioning**  
Relative positioning in model coordinates
- ✓ **Multiple model support**  
Up to 5 models can be displayed
- ✓ **Efficient workflow**  
Clone the last setup as starting point for a new measurement
- ✓ **Probe & data handling**  
Listing of active COM-ports, user-defined threshold on incoming data and addition of data to stored files

### Contact us for software updates or maintenance quotes:

- Continuous access to the latest updates and software features
- Expert technical support for installation, setup and trouble-shooting
- Multi-year discounts on request

Own an active license? Download the new software at [www.procap.tech/myprocap](http://www.procap.tech/myprocap)

streamwise gmbh  
Emil-Staub-Strasse 5  
8708 Männedorf | Switzerland

+41 44 545 32 50  
info@streamwise.ch